

Product Datasheet

**IKK- Alpha Antibody (Phospho-Ser176) : FITC
(orb2074041)**

Description

IKK- Alpha Antibody (Phospho-

Species/Host

Rabbit

Conjugation

FITC

Tested

ELISA, IHC, WB

Applications
Immunogen

The antiserum was produced against synthesized phospho-peptide derived from the Internal region of human IKK-alpha around the phosphorylation site of S176

Preservatives

Liquid. Purified antibody supplied in 1x PBS buffer.

Form/Appearance

Liquid. Purified antibody supplied in 1x PBS buffer.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: WB: 1:500 - 1:1000IHC: 1:50 - 1:100ELISA: 1:10000

Clonality

Polyclonal

MW

84 kDa

Uniprot ID
O15111, O14920
Expiration Date

12 months from date of receipt.