

## Product Datasheet

### Phospholipase A2 antibody (orb20745)

**Description**

Rabbit polyclonal antibody to Phospholipase

**Species/Host**

Rabbit

**Conjugation**

Unconjugated

**Preservatives**

Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran No preservative added, as it may interfere with the antibody activity. No foreign protein added.

**Form/Appearance**

Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran No preservative added, as it may interfere with the antibody activity. No foreign protein added.

**Storage**

Storage: The lyophilised PAb fraction is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product

**Note**

For research use only

**Clonality**

Polyclonal

**Source**

Phospholipase A2 isolated and purified from Crotalus durissus terrificus venom. Freund's complete adjuvant is used in the first step of the immunization procedure.

**Uniprot ID**

**P62022**

**Hazard Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in