

Product Datasheet

NMBR Antibody : HRP (orb2074658)

Description

NMBR Antibody : HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

ELISA, IF, WB

Applications
Immunogen

The antiserum was produced against synthesized peptide derived between 221-270 amino acids from the Internal region of human NMBR.

Preservatives

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Form/Appearance

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: WB: 1:500~1:1000; IF: 1:100~1:500; ELISA: 1:1000. Dilution Info: The optimal working dilution should be determined by the end user.

Clonality

Polyclonal

MW

43 kD

Uniprot ID
P28336
NCBI
NP_002502.2
Expiration Date

12 months from date of receipt.