

Product Datasheet

MRE11A Antibody : HRP (orb2076794)

Description

MRE11A Antibody : HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

ELISA, WB

Applications
Immunogen

The antiserum was produced against synthesized peptide derived between 431-480 amino acids from the Internal region of human MRE11A.

Preservatives

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Form/Appearance

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: WB: 1:500~1:3000; ELISA: 1:10000. Dilution Info: The optimal working dilution should be determined by the end user.

Clonality

Polyclonal

Uniprot ID
P49959
NCBI
NP_005581.2
Expiration Date

12 months from date of receipt.