

Product Datasheet

LATS2 Antibody : Biotin (orb2077321)

Description

LATS2 Antibody : Biotin

Species/Host

Rabbit

Conjugation

Biotin

Tested

ELISA, IF, IHC

Applications
Immunogen

The antiserum was produced against synthesized peptide derived between 541-590 amino acids from the Internal region of human LATS2.

Preservatives

Liquid. Purified antibody supplied in 1x PBS buffer.

Form/Appearance

Liquid. Purified antibody supplied in 1x PBS buffer.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: IHC: 1:50~1:100; IF: 1:100~1:500; ELISA: 1:20000. Dilution Info: The optimal working dilution should be determined by the end user.

Clonality

Polyclonal

Uniprot ID
[Q9NRM7](#)
NCBI
[NP_055387.2](#)
Expiration Date

12 months from date of receipt.