

Product Datasheet

THRA Antibody : Biotin (orb2078425)

Description

THRA Antibody : Biotin

Species/Host

Rabbit

Conjugation

Biotin

Tested Applications

ELISA, IHC, WB

Immunogen

The antiserum was produced against synthesized peptide derived between 10-59 amino acids from the N-Terminal region of human Thyroid Hormone Receptor alpha.

Preservatives

Liquid. Purified antibody supplied in 1x PBS buffer.

Form/Appearance

Liquid. Purified antibody supplied in 1x PBS buffer.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: WB: 1:500~1:3000; IHC: 1:50~1:100; ELISA: 1:10000. Dilution Info: The optimal working dilution should be determined by the end user.

Clonality

Polyclonal

Uniprot ID

P10827

NCBI

NP_001177847.1

Expiration Date

12 months from date of receipt.