

Product Datasheet

GAD2 Antibody : HRP (orb2078636)

Description

GAD2 Antibody : HRP

Species/Host

Rabbit

Conjugation

HRP

Tested Applications

ELISA, IHC, WB

Immunogen

The antiserum was produced against synthesized peptide derived between 545-594 amino acids from the C-Terminal region of human GAD1/2.

Preservatives

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Form/Appearance

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: WB: 1:500~1:3000; IHC: 1:50~1:100; ELISA: 1:5000. Dilution Info: The optimal working dilution should be determined by the end user.

Clonality

Polyclonal

Uniprot ID
[Q05329](#)
NCBI
[NP_001127838.1](#), [NP_000809.1](#)
Expiration Date

12 months from date of receipt.