



Product Datasheet

RNF93 antibody (orb2079)


Description	RNF93 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human RNF93/ZNF178 (101-200/465 aa)
Target	RNF93
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	Q9C019
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Stomach l...



Immunohistochemical staining of Rat kidn...



Immunohistochemical staining of rat test...