

Product Datasheet

NCF1 Antibody : Biotin (orb2079202)

Description

NCF1 Antibody : Biotin

Species/Host

Rabbit

Conjugation

Biotin

Tested

ELISA, IF, IHC

Applications

Immunogen

The antiserum was produced against synthesized peptide derived between 301-350 amino acids from the C-Terminal region of human p47 phox.

Preservatives

Liquid. Purified antibody supplied in 1x PBS buffer.

Form/Appearance

Liquid. Purified antibody supplied in 1x PBS buffer.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: IHC: 1:50~1:100; IF: 1:100~1:500; ELISA: 1:1000. Dilution Info: The optimal working dilution should be determined by the end user.

Clonality

Polyclonal

Uniprot ID

[P14598](#)

NCBI

[NP_000256.3](#)

Expiration Date

12 months from date of receipt.