

## Product Datasheet

### Pyruvate Kinase antibody (orb20824)

**Description**

Sheep polyclonal antibody to Pyruvate

<b>Species/Host</b>	Sheep
<b>Conjugation</b>	Unconjugated
<b>Preservatives</b>	Purified hyperimmune sheep IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added.
<b>Form/Appearance</b>	Purified hyperimmune sheep IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added.
<b>Storage</b>	Storage: The lyophilised PAb fraction is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Pyruvate kinase isolated and purified from rabbit muscle. Freund's complete adjuvant is used in the first step of the immunization procedure.
<b>Uniprot ID</b>	<b>P11974</b>
<b>Hazard Information</b>	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this