


Product Datasheet

Maltose Binding Protein antibody (orb2088)

Description	Maltose Binding Protein antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from E. coli Maltose Binding Protein (1-100/396 aa)
Target	MBP tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	P0AEX9
Dilution Range	WB=1:5000-20000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1:MBP-Tagged Fusion Pr...