

## Product Datasheet

### UBASH3A Antibody - middle region : HRP (orb2117636)

<b>Description</b>	UBASH3A Antibody - middle region : HRP
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	WB
<b>Immunogen</b>	The immunogen is a synthetic peptide directed towards the middle region of human UBASH3A
<b>Preservatives</b>	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
<b>Form/Appearance</b>	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
<b>Storage</b>	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Synthetic peptide located within the following region: PCSLPRRSRGIKDFENDPPLSSCGIFQSRIAGDALLDSGIRISSVFASPA
<b>Clonality</b>	Polyclonal
<b>MW</b>	74kDa
<b>Uniprot ID</b>	<b>P57075</b>
<b>NCBI</b>	<b>NP_061834</b>
<b>Expiration Date</b>	12 months from date of receipt.