


Product Datasheet

ALAS2 antibody (orb21226)

Description	ALAS2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human 5-aminolevulinate synthase, erythroid-specific, mitochondrial protein (50-587AA)
Target	Omega-gliadin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22557
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...



Immunohistochemical staining of human li...



Immunohistochemical staining of human th...