

Product Datasheet

C4 Binding Protein antibody (orb21238)

Description

Precipitating polyclonal sheep antiserum to

protein is a plasma...

Species/Host

Sheep

Conjugation

Unconjugated

Preservatives

Delipidated, heat inactivated, lyophilized, stable antiserum. No preservative added.

Form/Appearance

Delipidated, heat inactivated, lyophilized, stable antiserum. No preservative added.

Storage

Storage: The lyophilized antiserum is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Dilutions may be prepared by adding phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the antiserum. Diluted antiserum should be stored at +4°C, not refrozen, and preferably used the same day

Note

For research use only

Clonality

Polyclonal

Source

C4-binding protein is a plasma protein with a molecular weight of 500,000 and consists of 6 chains. Its concentration in serum is about 250 µg/ml. It is one of the inhibitors of the complement activation system. The protein combines up to six molecules of C4b at or near the C2 combining site and is thus preventing further association with C2. It has been isolated as a homogenous protein for use in the antiserum production. Freund's complete adjuvant is used in the first step of the immunization procedure.

Uniprot ID

P04003

Hazard Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.