



Product Datasheet

APEH antibody (orb21259)


Description	APEH antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Acylamino-acid-releasing enzyme protein (201-500AA)
Target	Guanine nucleotide-binding protein subunit beta-like protein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13798
Dilution Range	WB:1:2000-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human resected tissue



Immunohistochemical staining of human colon tissue



Immunofluorescence analysis of Hela cells