



Product Datasheet

ASAH2 antibody (orb21280)


Description	ASAH2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Neutral ceramidase protein (255-506AA)
Target	Pan Methyl Lysine
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NR71
Dilution Range	WB:1:500-5000, IHC-P:1:100-500, IF:1:50-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sk...



Western blot analysis of HeLa whole cell...



Immunohistochemical staining of human ki...