



Product Datasheet

ADRB2 (phospho-S346) antibody (orb213534)


Description	Rabbit polyclonal antibody to ADRB2 (phospho-S346)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S346 of human Beta-2 Adrenergic Receptor protein. The exact sequence is proprietary.
Target	ADRB2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P07550, P10608, P18762
Entrez	154, 11555
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Beta-2 Adrenergic...



Immunohistochemical analysis of Beta-2 A...



Immunofluorescent analysis of Beta-2 Adr...