



Product Datasheet

BNIP2 antibody (orb213611)

Description	Rabbit polyclonal antibody to BNIP2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human BNIP2. The exact sequence is proprietary.
Target	BNIP2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q12982 , O54940
Entrez	663 , 12175
Dilution Range	WB: 1-500-1-1000, IF/ICC: 1-100-1-500, FACS: 1-100-1-200
Expiration Date	12 months from date of receipt.



Western blot analysis of BNIP2 expression...



Immunofluorescent analysis of BNIP2 staining...