



Product Datasheet

CDK1/2/3 (phospho-T14) antibody (orb213710)


Description	Rabbit polyclonal antibody to CDK1
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T14 of human CDK1/2/3 protein. The exact sequence is proprietary.
Target	CDK1; CDK2; CDK3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P11440 , P97377 , Q00526 , Q80YP0 , Q63699 , P24941 , P06493 , P39951
Entrez	54237 , 12534 , 1017 , 12566
Dilution Range	WB: 1-500-1-1000, IF/ICC: 1-100-1-500
Expiration Date	12 months from date of receipt.



Western blot analysis of CDK1/2/3 (Phosp...



Immunohistochemical analysis of CDK1/2/3...



Immunofluorescent analysis of CDK1/2/3 (...)