



Product Datasheet

p27 Kip1 (phospho-T198) antibody (orb213715)

Description	Rabbit polyclonal antibody to CDKN1B
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T198 of human p27 Kip1 protein. The exact sequence is proprietary.
Target	CDKN1B
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P46414 , P46527
Entrez	1027 , 12576
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of p27 Kip1 (Phospho-T198) antibody.



Immunohistochemical analysis of p27 Kip1 (Phospho-T198) antibody.