


Product Datasheet

CHRNA7 antibody (orb213740)

Description	Rabbit polyclonal antibody to CHRNA7
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human CHRNA7. The exact sequence is proprietary.
Target	CHRNA7
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q494W8 , Q05941 , P36544 , P49582
Entrez	11441 , 1139 , 25302
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of CHRNA7 expression...