



Product Datasheet

CREB1 (phospho-S133) antibody (orb213775)


Description	Rabbit polyclonal antibody to CREB1 (phospho-S133)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S133 of human CREB protein. The exact sequence is proprietary.
Target	CREB1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P16220, Q01147, P15337
Entrez	1385, 12912, 81646
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200, IP: 1-10-1-100
Expiration Date	12 months from date of receipt.



Western blot analysis of CREB (Phospho-S133) antibody.



Immunohistochemical analysis of CREB (Phospho-S133) antibody.



Direct ELISA antibody dose-response curve.