



## Product Datasheet

### CRP antibody (orb213782)


|                            |                                                                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CRP                                                                                                                 |
| <b>Species/Host</b>        | Rabbit                                                                                                                                            |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                      |
| <b>Tested Applications</b> | IF, IH, WB                                                                                                                                        |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the center region of human C-Reactive Protein. The exact sequence is proprietary. |
| <b>Target</b>              | CRP                                                                                                                                               |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                         |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                         |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                                  |
| <b>Note</b>                | For research use only                                                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                                        |
| <b>Source</b>              | Rabbit                                                                                                                                            |
| <b>Uniprot ID</b>          | <a href="#">P02741</a> , <a href="#">P14847</a> , <a href="#">P48199</a>                                                                          |
| <b>Entrez</b>              | <a href="#">25419</a> , <a href="#">12944</a> , <a href="#">1401</a>                                                                              |
| <b>Dilution Range</b>      | WB: 1-500-1-1000                                                                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                   |



Western blot analysis of C-Reactive Prot...



Immunohistochemical analysis of C-Reacti...



Immunofluorescent analysis of C-Reactive...