



Product Datasheet

CTNNB1 antibody (orb213801)

Description	Rabbit polyclonal antibody to CTNNB1
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rabbit, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human Beta-catenin. The exact sequence is proprietary.
Target	CTNNB1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9WU82 , Q02248 , P35222
Entrez	84353 , 1499 , 12387
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Beta-catenin ex...



Immunohistochemical analysis of Beta-cat...