



Product Datasheet

E2F6 antibody (orb213868)


Description	Rabbit polyclonal antibody to E2F6
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IH, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human E2F6. The exact sequence is proprietary.
Target	E2F6
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O75461 , O54917
Entrez	1876 , 50496
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200, ChIP: 1-100-1-500
Expiration Date	12 months from date of receipt.



Western blot analysis of E2F6 expression...



Immunohistochemical analysis of E2F6 sta...



ChIP analysis of Mouse ES cell nuclear I...