



Product Datasheet

ETS1 (phospho-T38) antibody (orb213907)


Description	Rabbit polyclonal antibody to ETS1
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T38 of human ETS1 protein. The exact sequence is proprietary.
Target	ETS1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P14921 , P27577 , P41156
Entrez	23871 , 24356 , 2113
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200, IF/ICC: 1-100-1-500
Expiration Date	12 months from date of receipt.



Western blot analysis of ETS1 (Phospho-T38) antibody.



Immunohistochemical analysis of ETS1 (Phospho-T38) antibody.



Immunofluorescent analysis of ETS1 (Phospho-T38) antibody.