



Product Datasheet

FANCA (phospho-S1149) antibody (orb213919)

Description	Rabbit polyclonal antibody to FANCA
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S1149 of human FANCA protein. The exact sequence is proprietary.
Target	FANCA
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9JL70, O15360
Entrez	14087, 2175
Dilution Range	WB: 1-500-1-1000, IF/ICC: 1-100-1-500, IP: 1-10-1-100
Expiration Date	12 months from date of receipt.



Western blot analysis of FANCA (Phospho-S1149) antibody



Immunohistochemical analysis of FANCA (Phospho-S1149) antibody



Immunofluorescent analysis of FANCA (Phospho-S1149) antibody