



Product Datasheet

GRM2 antibody (orb214011)

Description	Rabbit polyclonal antibody to GRM2
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human mGLUR2. The exact sequence is proprietary.
Target	GRM2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q14416 , Q9QYS2 , P31422 , Q14BI2 , Q14832 , P31421
Entrez	2913 , 108068 , 2912 , 24415 , 108069 , 24416
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mGLUR2 expressi...



Immunohistochemical analysis of mGLUR2 S...