

## Product Datasheet

### IL1F5 ELISA Kit (Mouse) (orb2140211)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | IL1F5 ELISA Kit (Mouse)                                       |
| <b>Range</b>             | 0.031 - 2 ng/mL                                               |
| <b>Storage</b>           | Store as indicated in product manual.                         |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | Assay Info: Assay Methodology:<br>Quantitative Sandwich ELISA |
| <b>Sample Types</b>      | Serum, Plasma, Biological Fluids                              |
| <b>Uniprot ID</b>        | <b>Q9QYY1</b>                                                 |
| <b>Sensitivity</b>       | 18.75 pg/mL                                                   |
| <b>NCBI</b>              | <b>NP_001139559</b>                                           |
| <b>Expiration Date</b>   | Please enquire.                                               |