



Product Datasheet

HDAC1 (phospho-S421) antibody (orb214034)


Description	Rabbit polyclonal antibody to HDAC1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S421 of human Histone Deacetylase 1 protein. The exact sequence is proprietary.
Target	HDAC1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	009106, Q13547, Q4QQW4
Entrez	297893, 433759, 3065
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone Deacetylase 1 (HDAC1) protein.



Immunohistochemical analysis of Histone Deacetylase 1 (HDAC1) protein.



Immunofluorescent analysis of Histone Deacetylase 1 (HDAC1) protein.