



Product Datasheet

HNRNPK (phospho-S284) antibody (orb214045)

Description	Rabbit polyclonal antibody to HNRNPK
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S284 of human hnRNP K protein. The exact sequence is proprietary.
Target	HNRNPK
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P61979, P61980, P61978
Entrez	117282, 3190, 15387
Dilution Range	WB: 1-500-1-1000, IF/ICC: 1-100-1-500
Expiration Date	12 months from date of receipt.



Western blot analysis of hnRNP K (Phosph...



Immunohistochemical analysis of hnRNP K
...



Immunofluorescent analysis of hnRNP K (P...