


Product Datasheet


NRAS antibody (orb214056)

Description	Rabbit polyclonal antibody to NRAS
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human NRAS/HRAS/KRAS. The exact sequence is proprietary.
Target	NRAS; HRAS; KRAS
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P20171, P01116, Q04970, Q61411, P01112, P08556, P08644, P32883, P01111
Entrez	4893, 16653, 24605, 3845, 24525, 293621, 3265, 15461
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.




Western blot analysis of NRAS/HRAS/KRAS

...



Immunohistochemical analysis of NRAS/HRA...



Immunofluorescent analysis of NRAS/HRAS/...