



Product Datasheet

IGF1 R (phospho-Y1165/Y1166) antibody (orb214085)

Description	Rabbit polyclonal antibody to IGF1R
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y1165/Y1166 of human IGF1 Receptor protein. The exact sequence is proprietary.
Target	IGF1R
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P24062, P08069, Q60751
Entrez	25718, 16001, 3480
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of IGF1 Receptor (...)



Immunohistochemical analysis of IGF1 Rec...