



Product Datasheet

JUN antibody (orb214141)


Description	Rabbit polyclonal antibody to JUN
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, FC, IH, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human c-Jun. The exact sequence is proprietary.
Target	JUN
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P17325 , P05627 , P05412
Entrez	16476 , 24516 , 3725
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200, IP: 1-10-1-100, FACS: 1-100-1-300
Expiration Date	12 months from date of receipt.



Western blot analysis of c-Jun expression...



Immunohistochemical analysis of c-Jun staining...



Immunoprecipitation of c-Jun from 0.5 mg...