


Product Datasheet

SMAD3 (phospho-S204) antibody (orb214215)

Description	Rabbit polyclonal antibody against Smad3
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S204 of human SMAD3 protein. The exact sequence is proprietary.
Target	SMAD3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P84022, Q8BUN5, P84025
Entrez	4088, 25631, 17127
Dilution Range	WB: 1-500-1-1000, IF/ICC: 1-100-1-500, IP: 1-10-1-100
Expiration Date	12 months from date of receipt.




Western blot analysis of SMAD3 (Phospho-

...



Immunohistochemical analysis of SMAD3 (P...



Immunofluorescent analysis of SMAD3 (Pho...