



Product Datasheet

MYC (phospho-S62) antibody (orb214282)


Description	Rabbit polyclonal antibody to MYC
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S62 of human c-Myc protein. The exact sequence is proprietary.
Target	MYC
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P01108, P01106
Entrez	4609, 17869
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200, IF/ICC: 1-100-1-500
Expiration Date	12 months from date of receipt.



Western blot analysis of c-Myc (Phospho-S62) antibody.



Immunohistochemical analysis of c-Myc (Phospho-S62) antibody.



Immunofluorescent analysis of c-Myc (Phospho-S62) antibody.