



## Product Datasheet

### NFKB1 antibody (orb214308)


|                            |                                                                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to NFKB1                                                                                                               |
| <b>Species/Host</b>        | Rabbit                                                                                                                                            |
| <b>Reactivity</b>          | Human, Mouse, Porcine, Rat                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                      |
| <b>Tested Applications</b> | IF, IH, IP, WB                                                                                                                                    |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the center region of human NF-kappaB p105/p50. The exact sequence is proprietary. |
| <b>Target</b>              | NFKB1                                                                                                                                             |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                         |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                         |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                                  |
| <b>Note</b>                | For research use only                                                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                                        |
| <b>Source</b>              | Rabbit                                                                                                                                            |
| <b>Uniprot ID</b>          | <a href="#">Q63369</a> , <a href="#">P25799</a> , <a href="#">P19838</a>                                                                          |
| <b>Entrez</b>              | <a href="#">18033</a> , <a href="#">4790</a>                                                                                                      |
| <b>Dilution Range</b>      | WB: 1-500-1-1000, IHC-P: 1-100-1-200, IF/ICC: 1-100-1-500, IP: 1-10-1-100                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                   |



Western blot analysis of NF-kappaB p105/...



Immunohistochemical analysis of NF-kappa...



Immunofluorescent analysis of NF-kappaB ...