



Product Datasheet

PAK2 (phospho-S20) antibody (orb214356)

Description	Rabbit polyclonal antibody to PAK2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S20 of human PAK2 protein. The exact sequence is proprietary.
Target	PAK2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q64303 , Q8CIN4 , Q13177
Entrez	29432 , 224105 , 5062 , 100910732
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of PAK2 (Phospho-S20)



Immunohistochemical analysis of PAK2 (Phospho-S20)