



Product Datasheet

PFKFB2 (phospho-S483) antibody (orb214383)

Description	Rabbit polyclonal antibody to PFKFB2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S483 of human PFKFB2 protein. The exact sequence is proprietary.
Target	PFKFB2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P70265, Q9JJH5, O60825
Entrez	18640, 24640, 5208
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of PFKFB2 (Phospho-S483) antibody



Immunohistochemical analysis of PFKFB2 (Phospho-S483) antibody