



Product Datasheet

PLC gamma 2 (phospho-Y1217) antibody (orb214406)


Description	Rabbit polyclonal antibody to PLCG2
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y1217 of human PLC gamma 2 protein. The exact sequence is proprietary.
Target	PLCG2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P24135, P16885, Q8CIH5
Entrez	29337, 234779, 5336
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of PLC gamma 2 (Ph...



Immunohistochemical analysis of PLC gamm...



Immunofluorescent analysis of PLC gamma ...