

Product Datasheet

AICDA ELISA Kit (Human) (orb2144239)

Description	AICDA ELISA Kit (Human)
Range	78-5000pg/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue Homogenates
Uniprot ID	Q9GZX7
Sensitivity	39 pg/mL
NCBI	NP_001317272
Expiration Date	Please enquire.