



Product Datasheet

PRKACA antibody (orb214435)


Description	Rabbit polyclonal antibody to PRKACA
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human PKA C alpha/beta/gamma. The exact sequence is proprietary.
Target	PRKACA; PRKACB
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P68182 , P27791 , P22694 , P68181 , P05132 , P22612 , P17612
Entrez	293508 , 5567 , 5566 , 18747 , 18749
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200, IF/ICC: 1-100-1-500
Expiration Date	12 months from date of receipt.



Western blot analysis of PKA C alpha/beta/gamma.



Immunohistochemical analysis of PKA C alpha/beta/gamma.



Immunofluorescent analysis of PKA C alpha/beta/gamma.