


Product Datasheet

PKC zeta (phospho-T560) antibody (orb214448)

Description	Rabbit polyclonal antibody to PRKCZ
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T560 of human PKC zeta protein. The exact sequence is proprietary.
Target	PRKCZ
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q02956 , P09217 , Q05513
Entrez	5590 , 25522 , 18762
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of PKC zeta (Phospho-T560) antibody.



Immunohistochemical analysis of PKC zeta (Phospho-T560) antibody.