

Product Datasheet

JNK1/2/3 (phospho-T183) antibody (orb214454)

Description

Rabbit polyclonal antibody to MAPK8

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

IF, IH, WB

Applications

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T183 of human JNK1/2/3 protein. The exact sequence is proprietary.

Target

MAPK8; MAPK9; MAPK10

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

Q9WTU6, P49186, P49185, P45984, P49187, P53779, Q91Y86, P45983, Q61831

Entrez


116554, 25272, 26419, 26420, 5602, 5599, 50658, 5601

Dilution Range


WB: 1-500-1-1000, IHC-P: 1-100-1-200, IF/ICC: 1-100-1-500

Expiration Date


12 months from date of receipt.



Western blot analysis of JNK1/2/3 (Phosp...



Immunohistochemical analysis of JNK1/2/3...



Immunofluorescent analysis of JNK1/2/3 (...)