

Product Datasheet

CYP7B1 ELISA Kit (Human) (orb2144631)

Description	CYP7B1 ELISA Kit (Human)
Range	31.2-2000pg/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, plasma and other biological fluids
Uniprot ID	075881
Sensitivity	15.68 pg/mL
Expiration Date	Please enquire.