

Product Datasheet

PTGDS ELISA Kit (Mouse) (orb2145605)

Description	PTGDS ELISA Kit (Mouse)
Range	0.78-50ng/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, plasma and other biological fluids
Uniprot ID	O09114
Sensitivity	0.399 ng/mL
NCBI	NP_032989.2
Expiration Date	Please enquire.