



Product Datasheet

SOX9 (phospho-S181) antibody (orb214601)

Description	Rabbit polyclonal antibody to SOX9
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S181 of human SOX9 protein. The exact sequence is proprietary.
Target	SOX9
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q04887 , P48436
Entrez	20682 , 6662
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of SOX9 (Phospho-S181) antibody.



Immunohistochemical analysis of SOX9 (Phospho-S181) antibody.



Immunofluorescent analysis of SOX9 (Phospho-S181) antibody.