



Product Datasheet

THRA antibody (orb214669)

Description	Rabbit polyclonal antibody to THRA
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human THRA. The exact sequence is proprietary.
Target	THRA
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Based on IHC-P., the antibody THRA cross-react with THRA1, THRA2 and THRA3.
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P63058, P10827, P63059
Entrez	21833, 7067, 81812
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200
Expiration Date	12 months from date of receipt.



Western blot analysis of THRA expression...



Immunohistochemical analysis of THRA sta...