



Product Datasheet

p53 (phospho-S315) antibody (orb214692)


Description	Rabbit polyclonal antibody to p53 (phospho-S315)
Species/Host	Rabbit
Reactivity	Human, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S315 of human p53 protein. The exact sequence is proprietary.
Target	TP53
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P10361, P04637
Entrez	7157, 24842
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200, IF/ICC: 1-100-1-500, IP: 1-10-1-100
Expiration Date	12 months from date of receipt.



Western blot analysis of p53 (Phospho-S315) antibody.



Immunohistochemical analysis of p53 (Phospho-S315) antibody.



Immunofluorescent analysis of p53 (Phospho-S315) antibody.