



Product Datasheet

VEGFB antibody (orb214727)

Description	Rabbit polyclonal antibody to VEGFB
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human VEGFB. The exact sequence is proprietary.
Target	VEGFB
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P49765
Entrez	7423
Dilution Range	S-ELISA: 0.5 ug/ml, WB: 1-500-1-1000, IHC-P: 1-100-1-200
Expiration Date	12 months from date of receipt.



Western blot analysis of VEGFB expression...



Immunohistochemical analysis of VEGFB staining...